


# MIXPAC™ T-Mixer Colibri

## T-Mixer™ Colibri plus

for 18 to 75 ml Systems

- Impression of root canal
- impression taking of preparations and cavities
- reconstruction with root post & core build-up
- cementation of root post and post core




Ratio	Description	Outer Ø (mm)	Inner Ø (mm)	Gauge	Ele-ments	Ø (mm)	L (mm)*	Vol. (ml)	Order No.	Article
1:1, 2:1	mint, natural elements	1.1	0.8	19	12	4.2	63.8+15	0.62	160 701	MBT 4.2-12-N11

\* length mixing tip + length cannula

## T-Mixer™ Colibri

for 2.5 to 10 ml Systems

- Reconstruction with root post
- core build-up
- cementation of root post and post core
- bulk filling
- base/liner



Ratio	Description	Outer Ø (mm)	Inner Ø (mm)	Gauge	Ele-ments	Ø (mm)	L (mm)*	Vol. (ml)	Order No.	Article
1:1, 2:1	brown, natural elements	0.9	0.6	20	10	2.5	33.85+15	0.12	154 315	MLT 2.5-10-N09
1:1, 2:1	brown, natural elements	1.1	0.8	19	10	2.5	33.85+15	0.12	154 316	MLT 2.5-10-N11
1:1, 2:1	brown, natural elements	1.4	1.1	17	10	2.5	33.85+20	0.14	154 702	MLT 2.5-10-N14

\* length mixing tip + length cannula

\*\* Customer specific printing on packaging possible (Standard: Article-No., LOT-No., Made in Switzerland)

### Product advantages

### Customer benefits

Integrated cannula

No risk of cannula breaking off  
No additional intraoral tip required  
Saves time

360° rotatable cannula

Ergonomic access to areas which are difficult to reach

Bendable cannula with constant inner diameter

Constant flow of material even when the cannula is bent

Smoothly rounded medical grade stainless steel cannula

Safe application in the root canal and at the sulcus without risk of injury

